

Madison



A Dual motor Tilt in Space riser porter, with interchangeable pressure management seat and back cushions



Madison

The Madison, with dual tilt in space mechanism and rise function aims to encourage and support independence, but allow for ease of carer management.

The mechanism with a full range of reclined positions, as well as independent back recline and rise positions provides a wide range of positional options to support comfort and pressure management.

Key features at a glance:

- · Choice of interchangeable back cushions
- Choice of interchangeable pressure management seat cushions
- Dual Motor Tilt in Space rise and recline mechanism
- Seat depth adjustment with I" and 2" adjuster pads
- Heavy-duty push handle with rubber grips and a hook for the handset
- Seat and leg rest generously padded with memory foam and with a removable cover
- Simple to use handset with loop on arm
- · Four heavy duty lockable castors
- · Available in a wide range of fabrics.

Additional options:

- Sliding footplate
- Full length arm covers and replacement chaise cover and memory foam
- · Range of headrests
- Accupak rechargeable battery system for mains free operation.





"This chair is ideal for supporting those users who require encouragement to achieve a level of independence and mobility, whilst providing the adaptability to take account of changing conditions."

Kate SheehanIndependent Occupational Therapist









Tel: 0844 7766 001*

Email: info@reposefurniture.co.uk

Visit our web site for more details www.reposefurniture.co.uk